

Simple & Convenient Administration of Oxygen

Correct delivery of supplemental oxygen to a breathing victim will not cause harm in short - term use in most people (Australian Resuscitation Council, Guideline 10.4).

In nearly all injury situations, oxygen therapy is likely to have a favourable impact on the injury trauma.

OXI-port delivers oxygen at a constant rate of 8 ltrs/minute, which is the industry norm, so no need for first aiders to make decisions regarding what oxygen level to apply (as with Oxy-sok).

Oxygen devices should only be used by those who are trained in their use.

- ◇ Includes both a Shoulder-Strap and a Carry Handle.
- ◇ Does not include the oxygen cylinder.



OXI-port™
Leading the OXI-revolution

The latest and most user friendly oxygen equipment kit on the Australian market: the OXI-port™ 8000 leads the way. Features the latest KDK8000 OXI-flo™ regulator/flowmeter, 8L/Min fixed flow and incorporates Oxygen Flush. Designed to keep training to a minimum, but maximise the therapeutic effectiveness. Designed in Australia, for use by Australians.

Features: KDK OXI-flo™ regulator/flowmeter. Includes a carry bag with optional shoulder straps, Laerdal 'Pocket Mask' and instruction booklet. (Oxygen cylinder not included).

OxyPort™
SIMPLY ADD OXYGEN

© Medical Development Australia